

ZOLL[®]



GETTING TO
HIGH-QUALITY CPR

With Confidence, Consistency, and Results

Success Depends on Your CPR Quality

High-quality CPR is the cornerstone of a system of care that can optimize resuscitation outcomes. The 2013 CPR Quality Consensus Statement from the American Heart Association (AHA)¹ and the AHA and European Resuscitation Council 2010 Guidelines² underscore just how important CPR quality is. To deliver high-quality CPR, providers need to focus on depth, rate, limited interruptions, and release. ZOLL has solutions that are just right—enough information to improve CPR performance without being distracting—so that rescuers can reach and maintain proficiency in each of these four critical areas.

Depth

Guidelines recommend a depth of at least 2 inches, or 5.1 cm.

Compressions move blood through the body to protect vital organs. Adequate depth is required to essentially trap the heart between the sternum and spine and effectively squeeze the blood out.


How ZOLL technology helps: If a rescuer does not deliver compressions at an adequate depth, an audible and visual prompt to “Push Harder” will be initiated. Once proper depth is reached, it’s reinforced with a “Good Compressions” message. ZOLL’s professional defibrillators also display depth numerically.



Inter
Depth
Ra

Rate

Guidelines now recommend a compression rate between 100 and 120 compressions per minute. A rate of at least 100 compressions per minute is necessary to achieve perfusion sufficient to support minimal circulation, while a rate that is too fast often results in inadequate depth or leaning.




How ZOLL technology helps: An easy-to-follow metronome beep helps rescuers achieve the correct rate. ZOLL’s professional defibrillators also display the rate numerically.

Real CPR Help® technology is built into the ZOLL defibrillator electrodes, which include a single-use soft sensor that gathers CPR data and transfers it to the defibrillator. There are no extra pieces and nothing to remember to add. Just place the electrode and the CPR help is automatically enabled.

Interruptions

A key recommendation is to minimize interruptions in CPR, and the 2013 consensus statement urges rescuers to strive to increase the percentage of time compressions are delivered to 80%. Interruptions in chest compressions during CPR substantially reduce blood flow to the heart and brain. In addition, shock success is directly tied to pauses and begins to fall significantly after a pause as brief as 10 seconds.



How ZOLL technology helps: ZOLL's patented See-Thru CPR® technology filters the CPR artifact so rescuers can see if an organized rhythm develops, thereby minimizing the duration of pauses.

Interruptions
Release
te

Release

Guidelines emphasize that rescuers should not lean on the chest during compressions. The full upstroke, or release of a compression, is necessary to allow the heart to fill for the next compression. When rescuers fail to release the compression by not coming off the chest, pressure builds—making CPR ineffective.



How ZOLL technology helps: ZOLL provides a release indicator that shows the rescuer whether he or she is fully releasing and doing so fast enough to support filling the heart for the next down stroke. A "Release Fully" prompt reminds rescuers not to lean.

CPR Quality Review


All ZOLL products with Real CPR Help not only provide real-time CPR feedback but also record CPR performance data for easy download to ZOLL RescueNet® Code Review software for post-event review, analysis, and debriefing. Review is recommended as part of staff training and quality improvement and is considered in accreditation surveys from The Joint Commission. Code Review enables rescuers to see events as they occurred, providing ECG, vital signs, and the depth and rate of compressions for a full post-event analysis of the rescue. It automatically captures and organizes all resuscitation data for case-by-case review, quality assurance assessment, trending of program results, and attachment to patient electronic records. The comprehensiveness of this software is unmatched in the industry.



Get More from an Electrode

Rescuers shouldn't settle for "just electrodes" when they can use ZOLL's multifunction OneStep™, CPR Stat-padz®, and CPR-D-padz® electrodes.

The real advantage lies in the electrodes' CPR sensors, which enable real-time CPR feedback (Real CPR Help) without the need to add additional parts or pieces that are easily forgotten or left behind. The single-use sensor also reduces the possibility of transmitting infection from bodily fluids.



Code Review



- 1 Meaney PA, et al. *Circulation*. June 25, 2013; e-pub ahead of print.
- 2 2010 American Heart Association Guidelines for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care Science.

ZOLL Medical Corporation, an Asahi Kasei Group company, develops and markets medical devices and software solutions that help advance emergency care and save lives, while increasing clinical and operational efficiencies. With products for defibrillation and monitoring, circulation and CPR feedback, data management, fluid resuscitation, and therapeutic temperature management, ZOLL provides a comprehensive set of technologies that help clinicians, EMS and fire professionals, and lay rescuers treat victims needing resuscitation and critical care. For more information, visit www.zoll.com.

The Asahi Kasei Group is a diversified group of companies led by holding company Asahi Kasei Corp., with operations in the chemicals and fibers, homes and construction materials, electronics, and health care business sectors. Its health care operations include devices and systems for critical care, dialysis, therapeutic apheresis, transfusion, and manufacture of biotherapeutics, as well as pharmaceuticals, diagnostic reagents, and nutritional products. With more than 25,000 employees around the world, the Asahi Kasei Group serves customers in more than 100 countries. For more information, visit www.asahi-kasei.co.jp/asahi/en/.

ZOLL Medical Corporation Worldwide Headquarters

269 Mill Road
Chelmsford, MA 01824
978-421-9655
800-804-4356

For the ZOLL representative or distributor nearest you, visit www.zoll.com/contacts.